

#### Nathan A. Mall, MD, CIME

The Orthopedic Center of St. Louis Cartilage Restoration Center of St. Louis 14825 N. Outer 40, Suite 200; Chesterfield, MO 63017 P: 314.336.2555 F: 866.276.7014

# Post-Operative Protocol- MCL Repair/Reconstruction

Patient Name:	Date:	
Diagnosis:	Date of Surgery:	
Provide patient with home exercise program		

#### Week 1-2

- Ankle pumps every hour
- Post -op brace to maintain full extension.
- Quad sets & SLR (Brace on) with no lag
- TTWB with crutches
- Ice or Cryocuff Unit on knee for 20-30 minutes every hour
- Passive ROM exercises: Limits: 0 to 40 degrees.
- NO Hip adductor strengthening

### Week 3-4 (ROM 0-75deg, TTWB)

- Supervised PT 2- 3 times a week (may need to adjust based on insurance)
- Continue SLR's in brace with foot straight up, quad isometric sets, ankle pumps
- No weight bearing with knee in flexed position, TTWB with brace locked in full extension
- Patellar mobilization exercises
- Brace locked in full extension for ambulation and sleeping, and may unlock for sitting with limit 0-75deg.
- May **not** remove brace for HEP
- NO Hip adductor strengthening

### Week 5 (ROM as tolerated, TTWB)

- Continue with above exercises/ice treatments
- Advance ROM as tolerated with no limits with brace on
- Stationary bike for range of motion (short crank or high seat, no resistance) Ok to remove brace for bike here
- No weight bearing with knee in flexed position, continue TTWB with brace locked in full extension
- Perform scar message aggressively
- Progressive SLR program for quad strength with brace on start with 1 lb., progress 1 -2 lbs. per week
- Hamstring and hip PREs
- Seated leg extension (90 to 40 degrees) against gravity with no weight
- NO side lying Hip adductor strengthening

#### Week 6 (TTWB)

- Continue all exercises
- No weight bearing with knee in flexed position, TTWB with brace locked in full extension
- Flexion exercises seated AAROM
- AAROM (using good leg to assist) exercises (4-5x/day) with brace on
- Continue ROM stretching and overpressure into extension
- SLR's with brace on
- NO side lying Hip adductor strengthening
- Leg press 0-70 arc of motion

# The ORTHOPEDIC CENTER OF ST. LOUIS

#### Nathan A. Mall, MD, CIME

The Orthopedic Center of St. Louis Cartilage Restoration Center of St. Louis 14825 N. Outer 40, Suite 200; Chesterfield, MO 63017 P: 314.336.2555 F: 866.276.7014

### Week 7 (WBAT)

- Continue above exercises
- Start WBAT with brace on in full extension and D/C crutches when stable
- Hamstring and calf stretching
- Self ROM 4-5x/day using other leg to provide ROM
- Advance ROM as tolerated no limits, may remove brace for ROM
- Regular stationary bike if Flexion > 115
- · Heel raises with brace on
- Hip strengthening No side lying hip adduction

#### Week 8

- Continue above exercises
- Unlock brace for ambulation when quad control adequate
- Mini squats (0-60 degrees)
- 4 inch step ups
- Isotonic leg press (0 90 degrees)
- Lateral step out with theraband
- Hip strengthening

#### Week 9

• D/C brace if quad control adequate

Advance ROM, Goal: 0 to 115 degrees, walking with no limp

- Add ball squats
- Initiate retro treadmill with 3% incline (for quad control)
- Increase resistance on stationary bike
- Mini-squats and weight shifts
- Sport cord (bungee) walking
- 8 inch step ups
- 4 inch step downs

#### Week 10

- Begin resistance for open chain knee extension
- Swimming allowed, flutter kick only
- Bike outdoors, level surfaces only
- Progress balance and board throws
- Plyometric leg press
- 6-8 inch step downs
- Start slide board
- Jump down's (double stance landing)
- Progress to light running program and light sport specific drills if:

Quad strength > 75% contralateral side

Active ROM 0 to > 125 degrees

Functional hop test >70% contralateral side

Swelling < 1cm at joint line

No pain

Demonstrates good control on step down



# Nathan A. Mall, MD, CIME

The Orthopedic Center of St. Louis Cartilage Restoration Center of St. Louis 14825 N. Outer 40, Suite 200; Chesterfield, MO 63017 P: 314.336.2555 F: 866.276.7014

# Week 11-22

- Stairmaster machine
- If full ROM, quad strength > 80% contralateral side, functional hop test >85% contralateral side, satisfactory clinical exam:

Progress to home program for running. Progress to hops, jumps, cuts and sports specific drills. Begin to wean from supervised therapy.

4-5 months	
• Criteria to return to sports	
Full Active ROM	
Quadriceps >90% contralateral side	
Satisfactory clinical exam	
Functional hop test > 90% contralateral	
Other:	
Modalities:Electric StimulationUltrasoundHeat before/Ice after	
Frequency:x/ week x weeks	
	Nathan A Mall MI